





B U C C A L TUBES

RMO's NEW

FLI Series Buccal Tubes

deliver superior performance in an extremely small package. With Metal-Injection-Molded (MIM) construction for smooth comfortable contours, the ultra-low profile design features enhanced mesial openings that make wire insertion a snap. Combined with RMO's anatomical mesh bases, the FLI Series is ideal for both direct bond and molar band applications.





B U C C A L TUBES

Maxillary



FIRST MOLAR

Angulation	Torque	Distal Offset	M-D	R / L	.018"		.022"	
					On Bases	Tubes Only	On Bases	Tubes Only
0°	-10°	0°	4.5 mm	R	A08730	A06730	A08734	A06734
				L	A08731	A06731	A08735	A06735



SECOND MOLAR

Angulation	Torque	Distal Offset	M-D	R / L	.018"		.022"	
					On Bases	Tubes Only	On Bases	Tubes Only
0°	-10°	0°	3.5 mm	R	A08732	A06732	A08736	A06736
				L	A08733	A06733	A08737	A06737

Mandibular



FIRST MOLAR


Angulation	Torque	Distal Offset	M-D	R / L	.018"		.022"	
					On Bases	Tubes Only	On Bases	Tubes Only
0°	-25°	0°	4.5 mm	R	A08740	A06740	A08744	A06744
				L	A08741	A06741	A08745	A06745



SECOND MOLAR

Angulation	Torque	Distal Offset	M-D	R / L	.018"		.022"	
					On Bases	Tubes Only	On Bases	Tubes Only
0°	-25°	0°	3.5 mm	R	A08742	A06742	A08746	A06746
				L	A08743	A06743	A08747	A06747

ANATOMICAL MESH BASES



- FLI Series direct bond mesh bases incorporate specific contours and anatomical notches for maximum adaptation to tooth morphology
- Second molar bases feature distal-gingival relief for improved comfort on partially erupted molars



The World's Oldest
Synergistic, Bioprogressive®,
Breathing Enhancement
Orthodontic Company®

To order, please contact your RMO® Sales Representative or call 800.525.6375