


THE **FAST** TRACK TO STRAIGHT TEETH



AcceleDent
System

The AcceleDent System is a new device that enhances orthodontic force to safely accelerate tooth movement. The AcceleDent System is used for only 20 minutes each day, and will typically reduce treatment time by 30%. AcceleDent is complementary to all fixed appliances and aligners and uses a gentle pulsing force to increase the rate of alveolar bone remodeling.

AcceleDent is light, comfortable, hands-free, and can be used while multi-tasking or while engaged in a variety of other daily activities. Patients have reported use while in the shower,

Patient Benefits


- Reduced treatment time without compromised aesthetics
- Less likely to develop caries or gum disease with shortened treatment
- Clinical trial demonstrates excellent root resorption safety profile

Orthodontist Benefits

- Reduced treatment time enables improved practice productivity
- Innovative technology provides practice differentiation
- 100% compatible with any biomechanical system

AcceleDent System

- 1 Activator – small extraoral component that generates the vibrational force
- 2 Mouthpiece – connects to Activator; patient bites on to for transfer of pulsing force to the dentition
- 3 Charging Port – for recharging the device and providing usage compliance information
- 4 Travel Shell – allows the device to be easily and safely carried while “on the go”



To learn more about the clinically proven fast track to straight teeth, please visit www.aceledent.com or email info@orthoaccel.com



8275 El Rio
Suite 100
Houston, TX 77054
North America: 832.631.1660
Outside North America:
+44 (0) 20 3318 1915
www.aceledent.com

L05000 October 2009