



ENERGY CHAIN™



Endurance To Go The Distance

Energy Chain™ outperforms competitors start to finish. Its unique patented thermosetting resin formula provides superior elastic memory following initial stretch over competitor's plastic chain. University and clinical research demonstrate that RMO's Energy Chain delivers sustained light continuous force over a longer period than other available elastic chains. The tough patented material is stain-resistant, tasteless, non-toxic, and latex free.



Benefits:

- Patented formula provides superior shape memory following initial stretch
- Close and open space faster with consistent force up to eight weeks
- Reduce appointment intervals and save valuable chair time
- Less stress decay and less elongation over time compared to other available elastic chains
- Spacing variation between modules allows proper selection




Consider other products made of the same quality material as Energy Chain™:

Energy Thread™ Elastomeric Thread

RMO® Energy Thread™ is produced using the same material found in Energy Chain™. Energy Thread™ is extruded in two cross-section square sizes - .018" x .018" and .026" x .026." By producing Energy Thread™ in a square configuration, it ties snugly and eliminates knot slippage common with round thread.

- 20% more mass than round thread for enhanced force delivery
- Resists swelling, staining and discoloration
- Smooth surface and hydrophobic properties prevent plaque build-up better than round textured thread
- 25-feet per spool in clear or gray



Tough-O® Energy Ligatures

Made of the same great material as Energy Chain™, Tough-O® Energy Ligatures maintain a higher force between office visits than standard ligatures. Improved force longevity with Tough-O® ligatures allows for maximized treatment efficiency by improving slot prescription expression and decreasing arch wire slop.

- Same great material as Energy Chain™
- Decreased elasticity deformation
- High force retention

Be sure to pair your elastomeric products with RMO®'s exceptional line of brackets, bondable tubes, and wires.



Buying all your favorite RMO® products is just a fingertip away - visit our web store at www.rmortho.com

P00387 REV F

2165 Earlywood Drive, Franklin, IN 46131
 P 303.592.8200 F 303.592.8209 E rmosales@rmortho.com
 800.525.6375 | www.rmortho.com

